

We Set the Gold Standard for Breast Milk Shipping

Milk Stork is the first breast milk delivery service for business traveling moms, providing no-fuss express shipping or easy toting of refrigerated breast milk home to baby.

How it Works

Pump

According to your regular nursing schedule and refrigerate your milk.

Pack

Activate the cooling unit in your Milk Stork cooler. Then, pack your refrigerated milk in the breast milk storage bags provided.

Ship

All shipments come pre-labeled with a FedEx Priority Overnight label. Schedule a FedEx pickup on your own or ask your hotel to help.

Tote

Carry your cooler of milk home with you in the reusable tote provided.

For more information or to place your order, go to:
www.milkstork.com/Ohana



DOMESTIC

The Pump & Ship

Overnight your refrigerated breast milk home.

Available in Two Sizes:

34oz

For a mom who would like to overnight a daily supply of breast milk.

- 12" X 5" X 8"
- 72 hr minimum of refrigeration
- One day's worth of milk for most moms

72oz

For a mom who would like to overnight more milk less often.

- 12" X 8" X 8"
- 72 hr minimum of refrigeration
- Over two days' worth of milk for most moms

The Pump & Tote

Carry your refrigerated breast milk home.

Available in Two Sizes:

34oz

For a mom who would like to carry approximately a day's supply of breast milk home with her.

- Weighs 6.3 lbs when full
- 12" X 5" X 8"
- 60 hr maximum of refrigeration.
- One day's worth of milk for most moms

72oz

For a mom who would like to carry a multi-day supply of breast milk home with her.

- Weighs 9.6 lbs when full
- 12" X 8" X 8"
- 60 hr maximum of refrigeration
- Over two days' worth of milk for most moms

Disclaimers:

*FedEx Priority Overnight's shipping capabilities allow for next-day delivery to most destinations. FedEx does not operate on Sundays. Alaska, Hawaii and rural destinations may have longer transit times. Occasionally, weather and other events can cause disruptions to FedEx's service. Milk Stork cannot control for these disruptions. Milk Stork's Pump & Ship coolers provide at least 72hr of refrigeration.

**Whether carried on or checked as luggage, your Milk Stork cooler AND your breast milk will be inspected by TSA and may also be X-rayed. If you're toting your breast milk internationally, it is the user's responsibility to research and comply with each country's, airport's, and security agency's regulations for traveling with breast milk. Milk Stork is not responsible for inconsistent or incorrect enforcement of breast milk travel regulations and/or inspection of the Milk Stork coolers in the United States or internationally.



Questions?

support.milkstork.com

1-877-242-1306 | +1-510-356-0221

M-Sa, 6am - 8pm CT